

LabScat

Laboratory Turbidity Meter



Applications

- Laboratory turbidity measurement in bottles or cuvettes
- Verification of online measuring instruments
- Accelerated aging test in the bottle

Industries

- Beverage (beer, spirits, wine, soft drinks)
- Food
- Chemical industry
- Pharmaceutical industry

Advantages

- Compact design
- High measuring span
- Precise and reliable measurement also at 0 °C / +32 °F
- Colour-compensated dual-angle measurement according to MEBAK-Standard
- Bottle rotation and water bath to minimize interfering influences
- Monitoring of the water bath quality
- Integrated operation panel with colour touch screen display
- Quick verification with secondary standard
- Minimum maintenance, no tools required

